



PHAGO'SPORE

Phago'Spore is a biocidal disinfecting detergent for every surface including the surface of non-invasive medical devices.



Strengths

- Ready-for-use
- Bactericidal, fungicidal, virucidal and sporicidal (active in 10 minutes on C. difficile and in 5 minutes on Norovirus)
- Environmentally friendly (active substance peracetic acid)
- Food contact product (complies with the official Decision of September 6th 1999)

Packaging and Ordering

Unit	Units Per Carton	Product Number
750ml bottle	6	AIR_10-65003
5-litre canister	2	AIR_10-65007



Distributed Globally by Airande Pty Limited
 Level 13, Suite 1A, 465 Victori Avenue,
 CHATSWOOD NSW 2067 AUSTRALIA
 Tel: +61 2 9844 Fax: +61 2 9844
 E: sales@airande.global
 W: www.airande.global

Business Center - Cercle du Lac
 Rue de Rodeuhaie, 1
 1348 LOUVAIN-LA-NEUVE BELGIUM
 E: info.europe@airande.global



CHRISTEYS

Manufactured by Christeys France - Life Sciences
 31, Rue de la Maladrerie 44120 VERTOU France
 Tel: +33 (0) 2 40 80 27 27 Fax: +33 (0) 2 40 03 09 73
 W: www.phagogene.fr



PHAGO'SPORE

Phago'Spore is a biocidal disinfecting detergent for every surface including the surface of non-invasive medical devices.

APPLICATION

- PHAGO'SPORE is a biocide (TP02, TP04)
- PHAGO'SPORE is a disinfecting detergent for every surface including the surface of non-invasive medical devices.
- PHAGO'SPORE can be used on surfaces which are directly in contact with food if it is rinsed after application (complies with Decision of September, 8th 1999).
- PHAGO'SPORE is formulated from an active ingredient (peracetic acid) and guarantees a microbiological efficiency with a large spectrum, and ensuring a respect for the environment thanks to its biodegradability in acetic acid, oxygen and water.
- PHAGO'SPORE has a foam spray which avoids the formation of aerosols and provides hgreater protection to the end user.

COMPOSITION

- Active disinfectant:
 - Peracetic acid (CAS No. 79-21-0) 0.034% m/m
 - Hydrogen Peroxide (CAS No. 7722-84-1) 3.26% m/m
- Perfume
99.77% of the ingredients are natural in their origin

ORGANOLEPTIC & PHYSICAL-CHEMICAL CHARACTERISTICS

- Appearance: liquid
- Colour: colourless
- Perfume: Mint
- pH of pure product at 20°C: 2.00
- Density at 20°C: 1.000 - 1.020 g/cm²

HOW TO USE

For professional use only

- **PHAGO'SPORE** is ready for use. Do not dilute.
- PHAGO'SPORE can be used for disinfecting floors following the impregnation strips protocol (1L for 5 strips)
- Spray on a single-use linen and/or apply on the surface. Don't rinse.
- Don't rinse.*
- Let it work 5 minutes for a bactericidal activity or 15 minutes for a complete spectrum.
- *In case of disinfecting surfaces which may be in contact with food, rinse the surfaces with clean water (official Decree of September, 8th 1999)
- Use biocides carefully.
- Read the label and the product information before use.

RECOMMENDATIONS

- Use gloves and appropriate protective clothing.
- Do not let the product run off the surface being cleaned.
- Expiry date: 18 months from the manufacturing date.
- Consult the Safety Data Sheet on the following website:
www.quickfds.com

ANTIMICROBIAL ACTIVITY

Standard	Concentration	Time
BACTERICIDAL		
EN 1276 dirty conditions (<i>P. aeruginosa</i> , and <i>E. hirae</i>)	80%	1 min
EN 13697 dirty conditions (<i>P.aeruginosa</i> , <i>S.aureus</i> , <i>E.hirae</i> and <i>E.coli</i>)	10%	5 min
EN 13727 dirty conditions (<i>P.aeruginosa</i> , <i>S.aureus</i> and <i>E.hirae</i>)	8%	5 min
EN 16615 dirty conditions (<i>P. aeruginosa</i> , <i>S. aureus</i> , and <i>E. hirae</i>)	100%	5 min
FUNGICIDAL		
EN 13697 dirty conditions (<i>C.albicans</i> and <i>A. Niger</i>)	100%	15 min
EN 13624 dirty conditions (<i>C.albicans</i>)	50%	15 min
EN 13624 dirty conditions (<i>A. brasiliensis</i>)	80%	15 min
EN 1650 dirty conditions (<i>Aspergillus</i>)	80%	15 min
EN 16615 dirty conditions (<i>P. aeruginosa</i> , <i>S. aureus</i> , and <i>E. hirae</i>)	100%	5 min
VIRUCIDAL		
pREN 14476 : 2011 dirty conditions (<i>Poliovirus type 1</i>)	100%	15 min
pr EN 14476: 2011 dirty conditions (<i>Adenovirus type 5</i>)	100%	5 min
pr EN 14476: 2011 dirty conditions (<i>Norvirus murin</i>)	100%	5 min
SPORICIDAL		
EN 13704 clean conditions (<i>C. difficile</i>)	100%	10 min
EN 13704 clean conditions (<i>B. subtilis</i>)	80%	15 min
EN 13697 dirty conditions (<i>C. difficile</i> . sporulated form)	50%	15 min

Please refer to the Scientific File for more detailed information and reports.



Distributed Globally by Airandé Pty Limited
Level 13, Suite 1A, 465 Victori Avenue,
CHATSWOOD NSW 2067 AUSTRALIA
Tel: +61 2 9844 Fax: +61 2 9844
E: sales@airandé.global
W: www.airandé.global

Business Center - Cercle du Lac
Rue de Rodeuhaie, 1
1348 LOUVAIN-LA-NEUVE BELGIUM
E: info.europe@airandé.global



CHRISTEYNS

Manufactured by Christeyns France - Life Sciences
31. Rue de la Maladrerie 44120 VERTOU France
Tel: +33 (0) 2 40 80 27 27 Fax: +33 (0) 2 40 03 09 73
W: www.phagogene.fr

